

Technical Information Clorox Healthcare[®] Bleach Germicidal Wipes





OVERVIEW

Clorox Healthcare[®] Bleach Germicidal Wipes are engineered as a ready-to-use sporicidal cleaner-disinfectant system for healthcare facilities.

These EPA-registered disinfectant wipes contain sodium hypochlorite and other ingredients to kill *C. difficile* spores and other healthcare-relevant pathogens in 3 minutes or less on a variety of hard, nonporous healthcare surfaces.

INGREDIENT INFORMATION

Active ingredient: sodium hypochlorite.....0.55%* Other ingredients:.....99.45% Total (not including weight of wipe):.....100% *0.52% (5200ppm) available chlorine

AVAILABLE FORMATS

Clorox Healthcare[®] Bleach Germicidal Wipes are available in multiple sizes and formats to accomplish different disinfecting jobs.

At a Glance							
Clorox Healthcare [®] Bleach Germicidal Wipes							
EPA Registration Number	67619-12						
Product Classification	Pre-moistened wipes						
Dilution	Ready to use						
Disinfectant Type	1-step cleaner-disinfectant						
Active Ingredient Chemistry	Sodium Hypochlorite (Bleach) (1:10 dilution of bleach) 0.55% (5500 ppm)						
Odor	fruity, floral, bleach						
Shelf Life	Expiration date listed on each package: month, day, year						
Toxicity Signal Word	Caution						
HMIS Rating	1-0-1-A						
PPE Required	Follow your facility's protocols						

and formats to accomplish different disinfecting jobs.		Designed for	ltem #	Quantity	Shelf Life
Clinical Wipes (5"x6")		Daily Disinfection Small high-touch surfaces Medical equipment	30577	150 ct. canister, 6/case	12 months
Multi-Purpose Wipes (6.75"x9")		Daily Disinfection Patient room surfaces Medical equipment	35309	70 ct. canister, 6/case	14 months
			31469	20 ct. soft pack, 24/case	12 months
			31424	50 ct. carton, 6/case	12 months
Terminal Wipes (12"x12")	Discharge Disinfection	30358	110 ct. bucket, 2/case	12 months	
		Large areas	30359	110 ct. refill, 2/case	12 months



Technical Information Clorox Healthcare[®] Bleach Germicidal Wipes





DISINFECTION DIRECTIONS FOR USE[#]



 Don PPE, then REMOVE gross soil if visible.^{*} When disinfecting surfaces contaminated with *C. difficile* spores, always clean surfaces prior to disinfecting.



 WIPE surface with wipe until completely wet.



 To DISINFECT, allow surface to remain wet for contact time listed below: Bacteria - 30 seconds Viruses - 1 minute C. difficile spores - 3 minutes



 WIPE surface with a clean, damp cloth to prevent residue buildup.

Special Instructions for Cleaning Prior to Disinfection Against *Clostridium difficile* spores:

- Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks and eye covering.
- **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right on horizontal surfaces, and top to bottom on vertical surfaces to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.
- Infectious Materials Disposal: Materials used in the cleaning process that may contain feces/waste are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

Special Instructions for Using This Product to Clean and Decontaminate Against HIV-1, HBV and HCV on Surfaces/Objects Soiled with Blood/Body Fluids:

- **Personal Protection:** When handling items soiled with blood or body fluids, use disposable impervious gloves, gowns, masks and eye coverings.
- **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and other objects before applying this product.
- **Contact Time:** Allow surfaces to remain wet for 1 minute, let air dry. For all other organisms, see directions for contact times.
- **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi critical medical devices prior to sterilization or high-level disinfection.

** See product label for complete instructions.* This is a one-step product. Precleaning is required only if blood, vomit, fecal matter or other soil is visible, or in advance of disinfecting *C. difficile* spores.



Technical Information Clorox Healthcare® Bleach Germicidal Wipes



PRODUCTS & SAFETY

DISINFECTION EFFICACY CLAIMS

	Effective against	Kill Time
	Acinetobacter baumannii	30 sec
	Bordetella pertussis	30 sec
	Campylobacter jejuni	30 sec
	Clostridium difficile (spores)	3 min
	Enterobacter aerogenes	30 sec
	Enterobacter cloacae	30 sec
	Enterococcus faecalis (VRE) (Vancomycin-resistant)	30 sec
	Enterococcus faecalis (MDR) (Multidrug-resistant)	30 sec
	Enterococcus hirae	30 sec
	Escherichia coli	30 sec
	Escherichia coli 0157:H7	30 sec
	Escherichia coli ESBL	30 sec
	Escherichia coli (NDM-1)	30 sec
	Klebsiella pneumoniae	30 sec
	Klebsiella pneumoniae ESBL	30 sec
	Klebsiella pneumoniae (CRKP) (Carbapenem-resistant)	30 sec
	Klebsiella pneumoniae (NDM-1)	30 sec
bacteria	Legionella pneumophila	30 sec
ã	Listeria monocytogenes	30 sec
	Proteus mirabilis	30 sec
	Pseudomonas aeruginosa	30 sec
	Salmonella enterica	30 sec
	Serratia marcescens	30 sec
	Shigella dysenteriae	30 sec
	Staphylococcus aureus	30 sec
	Staphylococcus aureus (LRSA) (Linezolid-resistant)	30 sec
	Staphylococcus aureus (MRSA) (Methicillin-resistant)	30 sec
	Staphylococcus aureus (CA-MRSA) (Community-associated Methicillin-resistant)	30 sec
	Staphylococcus aureus (VRSA) (Vancomycin-resistant)	30 sec
	Staphylococcus aureus (VISA) (Vancomycin-intermediate-resistant)	30 sec
	Staphylococcus epidermidis (CoNS)	30 sec
	Streptococcus pneumoniae	30 sec
	Streptococcus pyogenes	30 sec
2	Mycobacterium tuberculosis (TB)	3 min

	Effective against	Kill Time			
	Adenovirus Type 2	1 min			
	Avian Influenza A Virus	1 min			
	Canine Parvovirus	3 min			
	Enterovirus D68	1 min			
	Feline Panleukopenia Virus	3 min			
	Hepatitis A Virus (HAV)	1 min			
	Hepatitis B Virus (HBV)*	1 min			
	Hepatitis C Virus* (HCV)	1 min			
	Herpes Simplex Virus (HSV-2)	1 min			
	Human Coronavirus	1 min			
3	Human Immunodeficiency Virus Type 1 (HIV-1)	30 sec			
A II US	Influenza A Virus	1 min			
	Influenza B Virus	1 min			
	Measles Virus	1 min			
	Middle East Respiratory Syndrome Coronavirus (MERS-CoV)	1 min			
	Norovirus*	1 min			
	Poliovirus Type 1 (Mahoney)	1 min			
	Respiratory Syncytial Virus (RSV)	1 min			
	Rhinovirus	1 min			
	Rotavirus	1 min			
	Severe Acute Respiratory Syndrome-Associated Coronavirus (SARS-CoV)	1 min			
19	Candida albicans	3 min			
	Candida glabrata	3 min			
	Trichophyton mentagrophytes (Athlete's Foot Fungus)	3 min			
	*Obtained via surrogate TOTAL: 58 Claims				

1. Tested in the presence of 3-part organic soil load following EPA Guidance for the Efficacy Evaluation of Products with Sporicidal Claims Against Clostridium difficile (June 2014)

Use products as directed. © 2017 Clorox Professional Products Company



Technical Information Clorox Healthcare[®] Bleach Germicidal <u>Wipes</u>



SAFETY INFORMATION

Reference the Safety Data Sheet (SDS) sheet for complete product safety information.

EPA Label Statements

This chemical is a pesticide product registered by the Environmental Protection Agency and is subject to certain labeling requirements under federal pesticide law. These requirements differ from the classification criteria and hazard information required for safety data sheets and for workplace labels of non-pesticide chemicals. Following is the hazard information as required on the pesticide label:

CAUTION: Liquid causes moderate eye irritation. Do not get in eyes or on clothing. Avoid contact with clothing. Wear protective eyewear. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. For sensitive skin or prolonged use, wear gloves.

FIRST AID: IF IN EYES: Hold eye open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment. QUESTIONS? Call (800) 234-7700.

INHALATION: Remove to fresh air. If breathing problems persist, see a doctor

PHYSICAL AND CHEMICAL HAZARDS: This product contains bleach. Do not use this product with other chemicals, such as ammonia, toilet bowl cleaners, rust removers or acid, as this releases hazardous gases.

Other Ratings

National Fire Protection Association (NFPA) Ratings

Health Hazard: Flammability: 0 Instability: 1 Physical and Chemical Hazards :none

Hazardous Materials Identification System (HMIS) Ratings

Health Hazard: 1 Flammability: 0 Physical Hazard: 1 Personal Protection: A