





OVERVIEW

Clorox Healthcare® Bleach Germicidal Cleaner is engineered as a ready-to-use sporicidal cleaner-disinfectant system for healthcare facilities.

This EPA-registered disinfectant contains sodium hypochlorite and other ingredients to kill *C. difficile* spores and other healthcare-relevant pathogens in 3 minutes or less on a variety of hard, nonporous healthcare surfaces.

INGREDIENT INFORMATION

Active ingredient: sodium hypochlorite	0.65%*
Other ingredients:	99.35%
Total:	100%
*0.639/ /617E nnm) available chloring	

^{*0.62% (6175} ppm) available chlorine

AVAILABLE FORMATS

Clorox Healthcare® Bleach Germicidal Cleaner is available in multiple sizes and formats to accomplish different disinfecting jobs.

At a Glance Clorox Healthcare® Bleach Germicidal Cleaner				
EPA Registration Number	56392-7			
Product Classification	Trigger spray			
Dilution	Ready-to-use			
Disinfectant Type	1-step cleaner-disinfectant			
Active Ingredient Chemistry	Sodium Hypochlorite (Bleach) (1:10 dilution of bleach) 0.65% (6500 ppm)			
Odor	fruity, floral, bleach			
Shelf Life	Expiration date listed on each package: month, day, year			
Toxicity Signal Word	Caution			
HMIS Rating	1-0-1-A			
PPE Required	Follow your facility's protocols			

disinfecting jobs.	·		Designed for	Item#	Quantity	Shelf Life
Trigger Sprays Convenient application Enable quick and easy surface coverage	R	7	• •	68967	22 oz. spray, 8/case	18 months
	68970	32 oz. spray, 6/case	18 months			
Pull-tops	3		Minimizing Aerosolization Leverage for larger surfaces or areas where sprays are not preferred	68832	32 oz. pull- top, 6/case	18 months
Refills		Easy refilling Minimize packaging waste and meet the needs of heavy users.	68973	64 oz. refill, 6/case	18 months	
Reillis			68978	128 oz. refill, 4/case	18 months	







DISINFECTION DIRECTIONS FOR USE#

TRIGGER SPRAY



 Don PPE, then REMOVE gross soil if visible. When disinfecting surfaces contaminated with *C. difficile* spores, always clean surface prior to disinfecting.*



2. SPRAY 6"-8" from surface until surface is thoroughly wet.



3. To DISINFECT, surface should remain wet for 1 MINUTE. To kill *C. difficile* spores, surface should remain wet for 3 MINUTES.



4. WIPE surface with a clean, damp cloth to prevent residue buildup.

Special Instructions for Cleaning Prior to Disinfection Against *Clostridium difficile* spores:

- Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks and eye covering.
- Cleaning Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right on horizontal surfaces, and top to bottom on vertical surfaces to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.
- Infectious Materials Disposal: Materials used in the cleaning process that may contain feces/waste are to be disposed of immediately in accordance with local regulations for infectious materials disposal.
- **Disinfection Procedure:** Spray 6"–8" from surface until surface is completely wet. Let stand for 3 minutes. Wipe with a clean, damp cloth or paper towel or allow to air dry.

Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV and HCV on surfaces/Objects Soiled with Blood/Body Fluids.

- **Personal Protection:** When handling items soiled with blood or body fluids, use disposable gloves, gowns, masks, and eye coverings.
- **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before applying this product.
- Contact Time: Let stand for 1 minute. For all other organisms, see directions for contact times.
- **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

^{*}See product label for complete instructions. This is a one-step product. Precleaning is required only if blood, vomit, fecal matter or other soil is visible, or in advance of disinfecting C. difficile spores.







DISINFECTION DIRECTIONS FOR USE#

PULL-TOP



 Don PPE, then REMOVE gross soil if visible. When disinfecting surfaces contaminated with *C. difficile* spores, always clean surface prior to disinfecting.*



2. APPLY directly onto hard, nonporous surface until surface is thoroughly wet.



3. To DISINFECT, surface should remain wet for 1 MINUTE. To kill *C. difficile* spores, surface should remain wet for 3 MINUTES.



4.WIPE surface with a clean, damp cloth to prevent residue buildup.

Special Instructions for Cleaning Prior to Disinfection Against *Clostridium difficile* spores:

- Personal Protection: Wear appropriate barrier protection such as gloves, gowns, masks and eye covering.
- Cleaning Procedure: Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with this product. Cleaning is to include vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right on horizontal surfaces, and top to bottom on vertical surfaces to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.
- Infectious Materials Disposal: Materials used in the cleaning process that may contain feces/waste are to be disposed of immediately in accordance with local regulations for infectious materials disposal.
- **Disinfection Procedure:** Spray 6"–8" from surface until surface is completely wet. Let stand for 3 minutes. Wipe with a clean, damp cloth or paper towel or allow to air dry.

Special Instructions for Cleaning and Decontamination Against HIV-1 (AIDS Virus), HBV and HCV on surfaces/Objects Soiled with Blood/Body Fluids.

- **Personal Protection:** When handling items soiled with blood or body fluids, use disposable gloves, gowns, masks, and eye coverings.
- **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before applying this product.
- Contact Time: Let stand for 1 minute. For all other organisms, see directions for contact times.
- **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

^{*}See product label for complete instructions. This is a one-step product. Precleaning is required only if blood, vomit, fecal matter or other soil is visible, or in advance of disinfecting C. difficile spores.







DISINFECTION EFFICACY CLAIMS

	Effective against	Kill Time
	Acinetobacter baumannii	1 min
	Bordetella pertussis	1 min
	Campylobacter jejuni	1 min
	Clostridium difficile (spores)	3 min
	Enterobacter aerogenes	1 min
	Enterobacter cloacae (NDM-1)	1 min
	Enterococcus faecalis (VRE) (Vancomycin-resistant)	1 min
	Enterococcus faecium (MDR) (Multidrug-resistant)	1 min
	Enterococcus faecium (VRE) (Vancomycin-resistant)	1 min
	Enterococcus hirae	1 min
	Escherichia coli ESBL	1 min
	Escherichia coli (NDM-1)	1 min
	Klebsiella pneumonia	1 min
o	Klebsiella pneumoniae ESBL	1 min
Bacteria	Klebsiella pneumoniae (Carbapenem-resistant)	1 min
Вас	Klebsiella pneumoniae (NDM-1)	1 min
	Proteus mirabilis	1 min
	Pseudomonas aeruginosa	1 min
	Salmonella enterica	1 min
	Serratia marcescens	1 min
	Staphylococcus aureus	1 min
	Staphylococcus aureus (MRSA) (Methicillin-resistant)	1 min
	Staphylococcus aureus (VRSA) (Vancomycin-resistant)	1 min
	Staphylococcus aureus (VISA) (Vancomycin intermediate resistant)	1 min
	Staphylococcus epidermidis (CoNS)	1 min
	Stenotrophomonas maltophilia	1 min
	Streptococcus pneumoniae	1 min
	Streptococcus pyogenes	1 min

	Effective against	Kill Time
	Adenovirus Type 2	1 min
	Avian Influenza A Virus	1 min
	Canine Parvovirus	1 min
	Enterovirus D68	1 min
	Feline Panleukopenia Virus	1 min
	Hepatitis A Virus (HAV)	1 min
	Hepatitis B Virus (HBV)*	1 min
	Hepatitis C Virus* (HCV)	1 min
	Herpes Simplex Virus (HSV-2)	1 min
	Human Immunodeficiency Virus Type 1 (HIV-1)	1 min
/iruses	Influenza A Virus	1 min
\$	Influenza A H1N1 Virus	1 min
	Influenza B Virus	1 min
	Measles Virus	1 min
	Middle East Respiratory Syndrome Coronavirus (MERS-CoV)	1 min
	Norovirus*	1 min
	Poliovirus Type 1 (Mahoney)	1 min
	Rhinovirus	1 min
	Rotavirus	1 min
	Severe acute respiratory syndrome (SARS)-associated coronavirus	1 min
1gi	Candida albicans	1 min
Ē	Trichophyton mentagrophytes (Athlete's Foot Fungus)	1 min

^{*}Obtained via surrogate

TOTAL: 50 Claims

Tested in the presence of 3-part organic soil load following EPA Guidance for the Efficacy Evaluation of Products with Sporicidal Claims Against Clostridium difficile (June 2014)







SAFETY INFORMATION

Reference the Safety Data Sheet (SDS) sheet for complete product safety information.

EPA Label Statements

This chemical is a pesticide product registered by the Environmental Protection Agency and is subject to certain labeling requirements under federal pesticide law. These requirements differ from the classification criteria and hazard information required for safety data sheets and for workplace labels of non-pesticide chemicals. Following is the hazard information as required on the pesticide label:

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling.

FIRST AID: If in Eyes: Hold eye open and rinse slowly and gently with water for 15–20 minutes. If present, remove contact lenses after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

PHYSICAL OR CHEMICAL HAZARDS: Do not use this product with ammonia or acids such as vinegar, rust removers, or toilet bowl cleaners. This product not intended for use on fabrics.

National Fire Protection Association (NFPA) Ratings

Health Hazard: 1 Flammability: 0 Instability: 1

Physical and Chemical Hazards: none

Hazardous Materials Identification System (HMIS) Ratings

Health Hazard: 1 Flammability: 0 Physical Hazard: 1

Personal Protection: A