

Technical Information

Clorox Healthcare® Fuzion® Cleaner Disinfectant



Products & Safety

Overview

Clorox Healthcare® Fuzion® Cleaner Disinfectant is engineered as a ready-to-use bactericidal, virucidal, fungicidal and sporicidal cleaner-disinfectant system for healthcare facilities.

This EPA-registered disinfectant contains sodium hypochlorite and other ingredients to kill *C. difficile* spores and 60 other pathogens in 2 minutes or less on a variety of hard, nonporous healthcare surfaces.

It leverages a revolutionary technology that delivers fast kill times with a neutralizer that breaks down the bleach to water and a small amount of salt to minimize residue and odor.

Available Formats

Clorox Healthcare® Fuzion® Cleaner Disinfectant is available as a trigger spray. It utilizes advanced nozzle technology to enable even surface coverage. It delivers larger droplets that reduce misting compared with a traditional swirl nozzle and will cling to vertical surfaces to reduce dripping.

At a Glance Clorox Healthcare® Fuzion® Cleaner Disinfectant


EPA Registration Number	67619-30
Product Classification	Trigger spray
Dilution	Ready-to-use
Disinfectant Type	1-step cleaner-disinfectant
Active Ingredient	Sodium Hypochlorite (Bleach) 0.39% (3900 ppm)
Odor	Fragranced
Expiration Date	Listed on each package: day, month, year
Toxicity Signal Word	Caution
HMIS Rating	0-0-0-X
PPE Required	Follow your facility's protocols

Ingredient Information

Active ingredient: sodium hypochlorite.....0.39%*

Other ingredients:.....99.61%

*yields 0.37% available chlorine

	Designed for	Item No.	Quantity	Shelf Life
Trigger Spray 	<ul style="list-style-type: none"> ▶ Facility-wide usage ▶ Enable quick and easy surface coverage 	31478	32 oz. spray, 9/case	24 months

Technical Information

Clorox Healthcare® Fuzion® Cleaner Disinfectant



Products & Safety



1. REMOVE gross soil if visible. For *C. difficile* spores or TB, always clean surface prior to disinfecting.*
2. (If trigger is in locked position) SQUEEZE the red tab under the nozzle and push downward to unlock the trigger.
3. Fully, rapidly depress the trigger lever until it physically touches the bottle 5-10 times until a full mist exits the nozzle.
4. SPRAY 6"-8" from surface until surface is completely wet.
5. To DISINFECT, let stand for 1 MINUTE. To kill *C. difficile* spores, allow 2 minutes contact. Then WIPE. Allow to air dry.

Special label instructions for cleaning prior to disinfection against *Clostridium difficile* spores:

- ▶ **Personal Protection:** Wear appropriate barrier protection such as gloves, gowns, masks and eye covering.
- ▶ **Cleaning Procedure:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop and/or sponge saturated with the sporicidal product. Cleaning includes vigorous wiping and/or scrubbing, until all visible soil is removed. Special attention is needed for high-touch surfaces. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading of the spores. Restrooms are to be cleaned last. Do not reuse soiled cloths.
- ▶ **Infectious Materials Disposal:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.
- ▶ **To Disinfect:** Spray 6-8 inches from surface until surface is thoroughly wet. Allow this product to remain wet for 2 minutes. Then wipe.

Special instructions for using this product to clean and decontaminate against HIV-1 and HBV on surfaces/objects soiled with blood/body fluids:

- ▶ **Personal Protection:** When handling items soiled with blood or body fluids, use disposable impervious gloves, gowns, masks and eye coverings.
- ▶ **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and other objects before applying this product.
- ▶ **Contact Time:** Allow surface to remain wet for HBV and HIV-1 for 1 min.
- ▶ **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to preclean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.

Technical Information

Clorox Healthcare® Fuzion® Cleaner Disinfectant



Products & Safety

Disinfection Efficacy Claims

Effective Against:
Spores - 2 min kill time
<i>Clostridioides difficile</i> spores
Bacteria - 1 min kill time
<i>Acinetobacter baumannii</i>
<i>Bordetella pertussis</i>
<i>Bordetella bronchiseptica</i>
<i>Burkholderia cepacia</i>
<i>Campylobacter jejuni</i>
<i>Enterobacter cloacae</i> (NDM-1)
<i>Enterococcus faecalis</i>
<i>Enterococcus faecalis</i> (VRE) (Vancomycin-resistant)
<i>Enterococcus faecium</i> (MDR) (Multidrug-resistant)
<i>Enterococcus hirae</i>
<i>Escherichia coli</i> O157:H7 (<i>E. coli</i>)
<i>Klebsiella aerogenes</i>
<i>Klebsiella pneumoniae</i> (CRE) (Carbapenem-resistant)
<i>Klebsiella pneumoniae</i> (NDM-1)
<i>Legionella pneumophila</i>
<i>Listeria monocytogenes</i>
<i>Proteus mirabilis</i>
<i>Pseudomonas aeruginosa</i>
<i>Salmonella enterica</i>
<i>Serratia marcescens</i>
<i>Shigella dysenteriae</i>
<i>Staphylococcus aureus</i>
<i>Staphylococcus aureus</i> (MRSA) (Methicillin-resistant)
<i>Staphylococcus aureus</i> (MDR) (Multidrug-resistant)
<i>Staphylococcus aureus</i> (VISA) (Vancomycin-intermediate)
<i>Stenotrophomonas maltophilia</i>
<i>Streptococcus pneumoniae</i>
<i>Streptococcus pneumoniae</i> (MDR) (Multidrug-resistant)
<i>Streptococcus pyogenes</i>

Viruses - 1 min kill time
Adenovirus Type 2
Adenovirus Type 5
Avian Influenza A Virus (H5N1)
Canine Adenovirus
Canine Parvovirus
Coxsackie Type B3 Virus
Enterovirus D68
Feline Panleukopenia Virus
Feline Rhinotracheitis Virus
Hepatitis A Virus
Hepatitis B Virus [†]
Hepatitis C Virus [*]
Herpes Simplex Virus Type 1 (HSV-1)
Human Coronavirus
Human Immunodeficiency Virus Type 1 (HIV-1)
Influenza A Virus
Influenza B Virus
Measles Virus
Middle East Respiratory Syndrome Coronavirus (MERS-CoV)
Murine Parvovirus
Norovirus [‡]
Poliovirus Type 1 (Chat)
Respiratory Syncytial Virus (RSV)
Rhinovirus
Rotavirus
SARS-Associated Coronavirus
SARS-CoV-2 (cause of COVID-19)
Fungi - 1 min kill time
<i>Candida albicans</i>
<i>Candida glabrata</i>
<i>Microsporum canis</i>
<i>Trichophyton interdigitale</i> (Athlete's Foot Fungus)
TB - 1 min kill time
<i>Mycobacterium bovis</i> (TB)
TOTAL: 66 CLAIMS
Use as directed on hard, nonporous surfaces.

* Obtained via Duck Hepatitis B Virus surrogate.

† Obtained via Bovine Viral Diarrhea Virus surrogate.

‡ Obtained via Feline Calicivirus surrogate.

Use product as directed. Clorox Healthcare recommendations. Your facility's protocol may vary.

Technical Information

Clorox Healthcare® Fuzion® Cleaner Disinfectant



Products & Safety

Safety Information

Reference the Safety Data Sheet (SDS) sheet for complete product safety information.

EPA Label Statements

This chemical is a pesticide product registered by the Environmental Protection Agency and is subject to certain labeling requirements under federal pesticide law. These requirements differ from the classification criteria and hazard information required for safety data sheets and for workplace labels of non-pesticide chemicals. Following is the hazard information as required on the pesticide label:

CAUTION: Moderate eye irritant. Do not get in eyes or on clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or going to the toilet.

FIRST AID: If in eyes: Hold eye open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for further treatment advice. Have the product container or label with you when calling a poison control center or doctor, or going for treatment.

Other Ratings

National Fire Protection Association (NFPA) Ratings

Health Hazard: 0

Flammability: 0

Instability: 0

Physical and Chemical Hazards: N/A

Hazardous Materials Identification System (HMIS) Ratings

Health Hazard: 0

Flammability: 0

Physical Hazard: 0

Personal Protection: X

Questions?

Call (800) 234-7700



Visit: www.CloroxPro.com

© 2023 Clorox Professional Products Company,
1221 Broadway, Oakland, CA 94612.
NI-60696