








Clorox Healthcare® Equipment Compatibility Guide

		Bleach					Hydrogen Peroxide	Quat
Furniture/Apparel								
		Clorox Healthcare® Bleach Germicidal Wipes	Clorox Healthcare® Bleach Germicidal Cleaners	Clorox Healthcare® Fuzion Cleaner Disinfectant	Dispatch® Hospital Cleaner Disinfectant Towels with Bleach	Clorox® Germicidal Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes	Clorox Healthcare® VersaSure Alcohol-Free Disinfectant Wipes
Manufacturer	Product							
Kwalu	All furniture products							
healthHcentric	IC+ Upholstery Solution™							
Prime Medical	SAF-T™ Scrubs							
Midmark	Upholstery							
	Manual Examination Tables							
DentalEZ Equipment Solutions	JV Generation Chair					1:10		
Graham-Field	LUMEX FR577RG/RGH Series Clinical Care Recliners					1:10		
Clinton Industries	Lab-X Series Padded Blood Drawing Chairs							
Activeaid, Inc.	Model 285/20 Tilt-In-Space Adult Shower/Commode/Wheel Chair							
Harloff Company, Inc.	Powder Coated and Stainless Steel Carts					0.05%		Clorox Disinfecting Wipes
Wheelchairs of Kansas	Model 38060000 Bariatric Lift and Transfer Unit							
A-Dec	Dental Chair Upholstery							

Legend

Product is mentioned by name in Cleaning & Care guide and is approved for use; if relevant, dilution is stated

Active ingredient/level is listed in Cleaning & Care guide as approved for use; if relevant, dilution is stated