

Clorox Healthcare Equipment Compability Guide



Ultrasound/Transducer Probes (Phillips)



Manufacturer	Product	Bleach				Hydrogen Peroxide	Quat Alcohol
		Clorox Healthcare® Bleach Germicidal Wipes	Clorox Healthcare® Bleach Germicidal Cleaners	Clorox Healthcare® Fuzion Cleaner Disinfectant	Dispatch® Hospital Cleaner Disinfectant Towels with Bleach	Clorox® Germicidal Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
Phillips	Ultrasound Transducers & Probes (System Names):						
	T = Transducer						
	C = Cable and Connector						
	H = Handle and Steering Mechanism						
	3D6-2 (HD11), 3D8-4, 3D9-3v (HD11, HD15, iU22), 7.5MHZ Endo	T/C	T/C	T/C	T/C	T/C	
	C5-1 (CX50, HD15, iE33, iU22), C5-2 (HD5, HD6, HD7, HD11, HD15, iE33, iU22), C6-2 Compact (CX30, Sparq), C6-3, C8-4v, C8-5, C8-5 Compact (CX30, CX50), C9-3io, C9-4, C9-4ec, C9-5ec	T/C	T/C	T/C	T/C	T/C	
	L8-4, L9-3, L11-3, L12-3 (HD5, HD6, HD7, HD11, HD15), L12-4 Compact (CX30), L12-5 50 Compact (CX50), L12 5 50 mm, L15-7io (HD11, HD15, iE33, iU22), L17-5	T/C	T/C	T/C	T/C	T/C	
	S3-1, S4-1 (iU22), S4-2, S4-2 Compact (CX30, Sparq), S5-1 (CX50, iE33, iU22), S5-2, S8, S8-3, S8-3 Compact (CX30, CX50, HD5), S12, S12-4, S12-4 Compact (CX50)	T/C	T/C	T/C	T/C	T/C	
	V6-2 (HD11XE, HD15, iU22), V8-4	T/C	T/C	T/C	T/C	T/C	
	X3-1, X5-1, X7-2 (iE33, iU22)	T/C	T/C	T/C	T/C	T/C	
	3D6-2 (iU22), 3D9-3v (Affiniti, ClearVue, EPIQ)	T/C	T/C	T/C	T/C	T/C	
	BP10-5ec	T/C	T/C	T/C	T/C	T/C	
	C5-1 (Affiniti, EPIQ), C5-2 - ClearVue, C6-2 Compact (Affiniti), C8-4v, C8-5 Compact (Affiniti, EPIQ), C9-2, C9-3v, C9-4v, C9-4v Compact, C10-3v, C10-4ec	T/C	T/C	T/C	T/C	T/C	
	L12-3 (Affiniti, CX50, EPIQ), L12-4 X7-2 (EPIQ), L12-4 Compact (Affiniti), L12 5 50 Compact (Affiniti, EPIQ), L15-7io (Affiniti, CX30, CX50, EPIQ), L15-7io (HD11, HD15, iE33, iU22), L18-5	T/C	T/C	T/C	T/C	T/C	
	S4-1 (ClearVue), S4-2 Compact (Affiniti), S5-1 (Affiniti, EPIQ), S8-3 Compact (Affiniti, EPIQ), S12-4 Compact (Affiniti, EPIQ)	T/C	T/C	T/C	T/C	T/C	
	V6-2 (Affiniti, ClearVue, EPIQ, HD5), VL13-5, VL13-5 Compact	T/C	T/C	T/C	T/C	T/C	
	X5-1 (EPIQ), X6-1	T/C	T/C	T/C	T/C	T/C	
	1.9MHz CW Pencil, 1.9MHz TCD	T/C	T/C	T/C	T/C	T/C	
	5.0MHz CW Pencil	T/C	T/C	T/C	T/C	T/C	
	D2cwc, D2tcd, D5cwc	T/C	T/C	T/C	T/C	T/C	
	TCD	T/C	T/C	T/C	T/C	T/C	
	S6-2mpt S7-3t, S7-2t Omni, S8-3t	C	C	C	C	T/C/H	
	OMNI III TEE	C	C	C	C	T/C/H	
	X7-2t	C	C	C	C	T/C/H	
	L10-4lap	T/C	T/C	T/C	T/C	T/C/H	
	15-6L	T/C	T/C	T/C	T/C	C	
	CO2 Transducer	M1460A					
		Avalon Fetal Transducer					

Legend	
	Product is mentioned by name in Cleaning & Care guide and is approved for use
	Active ingredient/level is listed in Cleaning & Care guide as approved for use, but product is not listed by name
	* Product mentioned in company's Cleaning & Care guide has the same EPA sub-registration