








Clorox Healthcare® Equipment Compatibility Guide

		Bleach					Hydrogen Peroxide	Quat
Ventilator Systems								
		Clorox Healthcare® Bleach Germicidal Wipes	Clorox Healthcare® Bleach Germicidal Cleaners	Clorox Healthcare® Fuzion Cleaner Disinfectant	Dispatch® Hospital Cleaner Disinfectant Towels with Bleach	Clorox® Germicidal Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes	Clorox Healthcare® VersaSure Alcohol-Free Disinfectant Wipes
Manufacturer	Product							
Covidien	Puritan Bennett							
	540 series (409 was OK'd for touchscreen and arm)							
	840 series (409 was OK'd for touchscreen and arm)							
	980 series (409 was OK'd for touchscreen and arm)							
Kendall SCD	700 Sequential Compression System							
AV Impluse	6060 Controller and Accessories							
Phillips	Triology 100 System							
	Respironics V60/Plus Ventilator					5%; 1:10		
Getinge-Maquet	SERVO i/s Ventilator System and Tube							

Legend

Product is mentioned by name in Cleaning & Care guide and is approved for use; if relevant, dilution is stated

Active ingredient/level is listed in Cleaning & Care guide as approved for use; if relevant, dilution is stated