



Frontline Staff Quick Reference Guide

*Disinfection Tool Kit
for the Operating Room
Using Hydrogen Peroxide Wipes*



The AORN Seal of Recognition was awarded to the Clorox Healthcare Disinfection Tool Kit for the Operating Room, Pre-Operative and PACU Areas on November 26, 2019 and does not imply that AORN approves or endorses any product or service mentioned in these materials. The AORN Seal of Recognition program does not include any CE credit for programs.





1. Stick to one cleaning path: clockwise or counterclockwise.
2. Use the Zone Approach to clean.
3. Clean from high to low.
4. Wipe from clean to dirty.

7 Steps for Success



PREPARE

Gather supplies, perform hand hygiene and put on PPE.



INSPECT

Inspect the surface to see if there is any visible soil or foreign material.



REMOVE SOIL (CLEAN)

If there is soil or foreign material on the surface, you will need to clean first with one wipe.



WIPE SURFACE (DISINFECT)

Using a fresh wipe, apply to the surface in an S-pattern, moving from clean to dirty and overlapping slightly on each sweep of the surface. Do not go over the same area twice with the same wipe. Use more than one wipe on large surfaces if necessary.



WAIT (CONTACT TIME)

Ensure the surface remains visibly wet for the full contact time listed on the product label.



DISCARD

Discard the used wipe when you have finished wiping the surfaces, or if the wipe becomes soiled or dry.



POLISH

Wipe surface with a clean, dry or damp cloth to prevent residue build-up and protect surfaces.

Hydrogen Peroxide Cleaner Disinfectant Wipes



Products & Safety

Product Overview

Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes

CONTACT TIME:*



DIRECTIONS FOR USE



1. Don PPE.
REMOVE gross,
visible soil.



2. APPLY directly
onto surface
evenly until
thoroughly wet.



3. To DISINFECT,
keep surface
wet for 1 minute.



4. POLISH
surfaces to
remove residue
if visible.

At a Glance Hydrogen Peroxide Cleaner Disinfectant Wipes

Active Ingredient	Hydrogen Peroxide
Format	Pre-moistened wipes
Dilution	Ready to use (no mixing required)
Contact time*	1 minute
Kills <i>C. difficile</i> spores?	No
Shelf Life	Expiration date listed on each package
PPE Required	Follow your facility's protocol

Refer to Safety Data Sheets in your Manager's Toolkit for complete information.

* Use as directed, refer to the EPA master label for complete list of organisms & contact times.

Use products as directed. Clorox Healthcare recommendations.
Your facility's protocol may vary.
© 2019 Clorox Professional Products Company.

NI-48555 Updated November 2019

Hydrogen Peroxide Cleaner Disinfectant Wipes



Products & Safety

Tips and Best Practices

About the product

- This is a one-step cleaner disinfectant product with hydrogen peroxide as the active ingredient.
- The formula also contains detergents and other ingredients to speed up the disinfecting power of hydrogen peroxide and help improve the product's cleaning power.
- These detergents are found in common cleaning products.
- While hydrogen peroxide itself has no smell, the other ingredients make the product smell like a common cleaner.




Residue and Surface Compatibility

- Hydrogen peroxide wipes are fast-acting and effective cleaners and disinfectants.
- Hydrogen peroxide disinfectants can be used facility-wide on a wide variety of surfaces including stainless steel and plastics.
- Remember that all disinfectants can leave some residue on surfaces.
- If you see some residue left behind, you can wipe it away with a clean, dry or damp cloth.

Operating Room: High-Touch Surfaces



Hydrogen Peroxide Wipes



CLEAN AND DISINFECT HIGH-TOUCH SURFACES WHICH INCLUDE:



Anesthesia cart



Overhead lights



Light switches



Vital signs monitors



Mayo stands



**Door handles/
push plates**



Infusion pump



Supply carts



Anesthesia Machine



Computer on wheels



IV pole



Surgical stool



**OR table, straps,
positioning/
transfer devices**



Back table




**Trash can/
kick bucket**

Operating Room Cleaning & Disinfecting: Before First Case



Protocols

Hydrogen Peroxide Wipes




BEFORE FIRST CASE / START OF DAY CLEANING

- Nursing staff should damp dust all horizontal/level surfaces and equipment in the OR using a clean, lint-free or microfiber cloth moistened with an EPA-registered healthcare disinfectant.
- Equipment brought from outside of the OR must be cleaned and disinfected before being brought into the OR.

Operating Room Cleaning & Disinfecting: Between Cases



Hydrogen Peroxide Wipes



BETWEEN CASES

- Wear the appropriate PPE including gloves per facility protocol.
- Clean to remove gross soil before disinfecting.
- Disinfect equipment in the room. Clean from high to low, and from clean to dirty.
- Use a new wipe for each piece of equipment. You may need to use more than one wipe on large equipment such as the OR table.
- Keep the surface visibly wet for the full contact time.
- Clean and disinfect the floor in the immediate patient area if soiled.
- Allow all surfaces to dry before they come in contact with the next patient.

Operating Room Between-Case Cleaning Checklist




Checklist

After Every Case		
Cleaned	Not Applicable	Surface/ Piece of Equipment
		Anesthesia machine, cart, IV pump & pole
		Back table
		Chairs, stools, step stools
		Computer (including keyboard, mouse, screens)
		Electronic communication devices (e.g., cell phone, pager)
		Equipment & booms
		Floors
		Mayo & ring stands
		OR light & handle
		OR table including controls, mattress, side rails, table straps, attachments
		Positioning and transfer devices
		Patient monitors
		Regulators (e.g., suction, medical gases)
		Waste receptacles & kick buckets
		High-touch surfaces (e.g., door handle or push plates, light switches, telephone, etc.)
Only If Items Are Soiled		
		Cabinet doors
		Linen hamper
		Radiology equipment
		Scrub sinks
		Supply carts
		Vents/grills for return air
		Walls & ceilings
		Wheels & casters

Operating Room Cleaning & Disinfecting: Terminal



Hydrogen Peroxide Wipes



END OF DAY / TERMINAL CLEANING

- Start cleaning after the last surgical case of the day.
- Follow the procedures for cleaning between cases.
 - See Responsibility Grid for additional items/surfaces to be cleaned and disinfected.
- Move standing equipment to clean and disinfect the entire operating room floor.
- Clean and disinfect the walls and ceiling if necessary; clean and disinfect the floor per protocol.
- Clean casters and wheels on portable equipment.
- Disassemble, clean and disinfect, and dry all cleaning equipment before reuse and storage.

Operating Room Terminal Cleaning Checklist



Each Day Room Is Used

Cleaned	Not Applicable	Item	Comments
		Anesthesia machine, cart, IV pump & pole	
		Back table	
		Cabinet doors	
		Chairs, stools, step stools	
		Computer (inc. keyboard, mouse, screen)	
		Equipment & booms	
		Floors	
		High-touch surfaces (e.g., door handle or push plates, light switches, telephone)	
		Linen hamper	
		Mayo & ring stands	
		OR light & handle	
		OR table including controls, mattress, side rails, table straps, positioning and transfer devices, attachments	
		Patient monitors	
		Positioning and transfer devices	
		Radiology equipment	
		Regulators (e.g., suction, medical gasses)	
		Scrub sinks	
		Sterile processing areas & department	
		Supply carts	
		Vents/grills for return air	
		Walls & ceilings	
		Waste receptacles & kick buckets	
		Wheels & casters	