

Between-Procedure Cleaning In the OPERATING ROOM



Product Information


Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes



1 minute contact time*




* Use as directed, refer to the EPA master label for complete list of organisms & contact times.




1. Don PPE. **REMOVE** gross, visible soil.
2. **WIPE** surface with wipe until thoroughly wet.
3. To **DISINFECT**, keep surface visibly wet for wet for 1 minute.
4. **POLISH** surfaces to remove residue if visible


Between-Procedure Cleaning Protocol



Cleaning Sequence



Key Surfaces for Disinfection



The AORN Seal of Recognition was awarded to the *Clorox Healthcare Disinfection Tool Kit for the Operating Room, Pre-Operative and PACU Areas* on November 26, 2019 and does not imply that AORN approves or endorses any product or service mentioned in these materials. The AORN Seal of Recognition program does not include any CE credit for programs.